

## Raz-SP Self Propel Rehab Shower Commode

The RAZ-SP (Self-Propel) rehab shower commode chair comes standard with 4" (10 cm) of height adjustment, tension-adjustable back upholstery, self-propelling rear wheels, 5" (125 mm) duallocking casters, flip-up padded armrests, and commode pan. The easy-to-clean, stainless steel frame is designed with unobstructed side and front access for the user or attendant.



## **FEATURES**



SPECIFICATIONS	Raz-SP 22" wheels	Raz-SP 24" wheels
WEIGHT CAPACITY	350 lbs (160 kg)	350 lbs (160 kg)
FRAME WIDTH (seat and back frame)	Standard: 18" (46 cm)	330 (33 (33 ),
OVERALL WIDTH (with standard armpad)	26.5" (67 cm)	26.5" (67 cm)
OVERALL LENGTH (with footrest)	39" (99 cm)	39" (99 cm)
SEAT WIDTH (commode seat)	18" (46 cm)	18" (46 cm)
SEAT DEPTH (front of back post to front of seat)	16" or 19" (40.5 or 48 cm)	16" or 19" (40.5 or 48 cm)
CLEARANCE HEIGHT UNDER SEAT	16", 17", 18", 20" (41 cm, 43 cm, 46 cm, 51 cm) without pan 15.25", 16.25", 17.25", 19.25" (39 cm, 41 cm, 44 cm, 49 cm) with pan	15.75", 16.75", 17.75", 19.75" (40 cm, 43 cm, 45 cm, 50 cm) without pan brackets 15", 16", 17", 19" (38 cm, 41 cm, 43 cm, 48 cm) with pan brackets
ARMREST HEIGHT ABOVE SEAT (17" Adjusta-Back)	7" - 10" in $\frac{1}{2}$ " increments (18 - 25 cm)	7" - 10" in ½" increments (18 - 25 cm)
SEAT-TO-FLOOR HEIGHT	20" - 24" in 1" increments (51 - 61 cm)	20" - 24" in 1" increments (51 - 61 cm)
PRODUCT WEIGHT (with footrests)	47 lbs (21.3 kg)	49 lbs (22kg)
FRAME AND HARDWARE MATERIAL	Stainless Steel (Lifetime warranty against corrosion)	Stainless Steel (Lifetime warranty against corrosion)

All measurements taken with standard molded seat and standard arm pad.

Two-year limited warranty covers the frame against manufacturing defects. A one-year limited warranty covers all other components, subject to normal wear and tear exclusions.

Specifications, and options are subject to change without prior notice.

## **MANUFACTURED BY:**

RAZ Design Inc. 22 Howden Road, Toronto, Ontario Canada, M1R 3C7 T: (877) 720-5678 E: info@razdesigninc.com

Raz Design products are distributed in Australia by:

**Paragon Mobility** 

3/10 Tapnor Crescent, Brendale, QLD, 4500 Australia Tel: 13000 652 382

Fax: (07) 3112 6125 www.paragonmobility.com.au

