







The K Care Electric Bath Lift is designed to assist users to safely transfer into and out of a bath.

Electric Bath Lift KA610

Using a bath lift can enhance safety and independence in the bathroom, allowing individuals to enjoy a bath whilst minimising assistance from others.

Specifications

KA610

1071mm Overall Height Seat Height Range 60-480mm Overall Width 360mm Seat Width 350mm 430mm Seat Depth Overall Length - Folded 728mm Overall Length - Reclined 950mm 140° Backrest Angle Weight Without Accessories 11kg Max User Weight 140kg

Important: If you or someone you know is considering a bath lift, it's advisable to consult with a healthcare professional or occupational therapist to determine the most suitable option based on individual requirements.



Features

- · Easy to assemble and use.
- Compact in design with a folding frame for easy storage and transportation.
- High-quality durable seating pad and backrest.
- Easily positioned in standard-sized bathtubs (non-standard or specially customised bathtub sizes may not be suitable).
- · Suction cups for secure positioning.
- The bath lift can be raised and lowered during bathing.
- Fold out side flaps assist with transfer to and from the bathlift.
- The backrest can be locked to prevent accidental tilting.

- The bath lift is waterproof (except for the charger).
- The bath lift is controlled with a hand set

 using a rechargeable battery that is
 integrated into the handset. The handset
 includes a battery indicator, up/down
 controls and an emergency stop button.
- Australian Register of Therapeutic Goods – ARTG no 390530.
- Independently tested to International standards for Product Safety, Electrical Safety, and EMC.
- Power adaptor has Australian Electrical Compliance.





